

Technical Datasheet

Product Information

Article Description	:	WMC non-woven tissue rolls Medi Wipes, 90 sheets
Article Number	:	3050716 C
EAN/GTIN	:	4260271586290
Material	:	hydroentangled fleece, 100% PET, 60 g/m ²
Units / Carton	:	6 rolls à 90 sheets unfolded with pull-out function
Tissue Dimensions	:	± 200 x 360 mm
Subpackaging	:	individually shrink-wrapped in foil

Logistic Information

Carton Measures	:	± 400 x 270 x 225 mm
Cartons / Pallet	:	63 cartons
Layers / Pallet	:	7 layers with 9 cartons each
Carton Weight	:	± 2,69 kg
Pallet Measures	:	± 1200 x 800 x 2000 mm (incl. pallet)
Pallet Weight	:	± 188 kg (incl. pallet)



The data only describes product-relevant information and is based on our current state of knowledge. They do not constitute a guarantee of product properties within the meaning of the statutory warranty regulations.

Stand: 26.05.2025 | STV