

**AGANI™**

*Hypodermic needle*



# Product Solutions

**Your needles of choice**  
for a wide range of applications



 **TERUMO**

The safety and  
comfort of your patients  
is in your hands.




**Injection & Access**



# Smooth and easy puncture increases patient comfort.

Make the administration of medication as painless as possible with the special geometry of the sharp lancet point and silicone-coated cannula that allows a smooth penetration and prevents any knocking effect.<sup>1</sup>



Precisely shaped needle point  
with silicone coating.

Agani hypodermic needles have a sharp lancet tip and an outside layer of silicone lubricant.

**The sharp  
point**

- Facilitates the penetration of the skin
- Helps to minimize the trauma

**The silicone  
coating**

- Provides good gliding properties
- Allows a smooth insertion

# Secure and easy handling

increases safety and user comfort.



The steeper point angle of the short bevel<sup>2</sup> shortens the length of the needle bevel, thereby reducing the risk of vascular perforation.

The outside diameter of the needle can be recognized by the ISO colour-coded hub.

- Easy identification of the required size
- Easy recognition of flashback through the semi-transparent hub

All needles come sterile, individually packed in color-coded blister packaging.

- Supplied in boxes of 100 needles
- Comprehensive information on the packaging



# Choose from a wide range from 18G up to 30G.



## 18G Cannula (ø 1.2 mm)

Length 38 mm  
Colour Pink



## 19G Cannula (ø 1.1 mm)

Length 25, 38, 50 mm  
Colour Cream



## 20G Cannula (ø 0.9 mm)

Length 25, 38 mm  
Colour Yellow



## 21G Cannula (ø 0.8 mm)

Length 16, 25, 38, 50 mm  
Colour Green



## 22G Cannula (ø 0.7mm)

Length 32, 38 mm  
Colour Black



## 23G Cannula (ø 0.6 mm)

Length 16, 25, 32 mm  
Colour Blue



## 24G Cannula (ø 0.55 mm)

Length 25 mm  
Colour Purple



## 25G Cannula (ø 0.5 mm)

Length 16, 25 mm  
Colour Orange



## 26G Cannula (ø 0.45 mm)

Length 13, 23 mm  
Colour Brown




## 27G Cannula (ø 0.4 mm)

Length 16, 19 mm  
Colour Grey



## 30G Cannula (ø 0.3 mm)

Length 13 mm  
Colour Yellow

CE 0123 Does not contain components made of PVC or natural rubber latex. 

For further ordering information, please refer to the AGANI hypodermic needle product information card.

# AGANI™

Hypodermic needle



## Life is easier with needles from Terumo.

Agani hypodermic needles provide you and your patients with the best possible levels of care and comfort.

Sharp lancet point – easy puncture

Silicone coating – smooth needle movement

Thin walled needles – increased flow rate

Short bevel options – reduced risk of vascular perforation

Large portfolio – suitable for a wide range of applications

### Peak performance even in little things – made easy!



**Terumo Corporation**  
+81 3 3374 8111

**Terumo Europe NV**  
+32 16 38 12 11

Your Terumo partner

[www.terumo-europe.com](http://www.terumo-europe.com)

© Registered Trademark  
Published by Terumo Europe NV

MP210GB-0919FK-1(10-19)E